

MODEL 580



- | | |
|------|---|
| ● UX | Ceramic microphone for moderate losses. |
| ● S | For severe hearing losses. |
| ● U | For moderate losses. |
| ● PC | For minor losses and nerve cases. |

HEARING FOR ACTIVE PEOPLE

Model 580 means hearing for active people, because of its small size and long-life circuit. This slim behind-the-ear aid is contoured to fit behind the smallest ear. It is easy to operate with combined switch and volume control. The new transistor technique has also been incorporated in this series and makes it possible to achieve up to 500 hours of listening with the small standard battery 675 or equivalent. Accumulator cell DK 20 can also be used.

The series is available in four versions, all in the same small size, and all with Oticon's excellent tonal quality.